


Experimental Planning



Experimental planning considerations

1. Replicates
2. Confounding factors
3. Batch effects


Experimental planning considerations: Replicates

1. Replicates

- **Technical replicates**: Illumina has low technical variation unlike microarrays, hence technical replicates are unnecessary.
- **Biological replicates**, are absolutely essential. Have at least 3!
 - For **cultured cells** try to incorporate as much variation as possible
 - Avoid **pooling** RNA from multiple biological replicates

Experimental planning considerations: Replicates

More replicates yield more accurate analyses:



Liu, &., et.al., *Bioinformatics* (2014) 30(3): 301-304, <https://doi.org/10.1093/bioinformatics/btt688>

Experimental planning considerations: Replicates

Replicates versus sequencing depth:

- **Replicates** are almost always preferred to greater sequencing depth for bulk RNA-Seq.
- General **differential gene expression**:
 - ENCODE guidelines suggest **30 million SE reads** per sample (stranded).
 - **15 million reads** per sample is often sufficient, if adequate number of replicates.
 - Spend money on more biological replicates, if possible.

Experimental planning considerations: Replicates

Replicates versus sequencing depth:

- Detection of **lowly-expressed genes**:
 - Similarly benefits from replicates more than sequencing depth.
 - However, would likely want to sequence deeper with at least 30-60 million reads.

Experimental planning considerations: Replicates

Replicates versus sequencing depth:

- Isoform-level differential expression:
 - More depth (> 60 million reads per sample) and paired-end reads are suggested.
 - However, choose biological replicates over paired/deeper sequencing.
 - Perform careful QC of RNA quality. Be careful to use high quality preparation methods and restrict analysis to high quality RIN # samples.



Experimental planning considerations: Replicates

Replicates versus sequencing depth:

- Other types of RNA analyses (intron retention, small RNA-Seq, etc.) will have different suggestions.
 - But almost always more biological replicates is better!

Experimental planning considerations: Confounding


2. Confounding factors



Experimental planning considerations: Confounding

Avoiding confounding:



- Ensure animals in each condition are all the **same sex, age, litter, and batch**, if possible
- If not possible, then ensure to split the animals equally between conditions



Experimental planning considerations: Batches

3. Batch effects

Batch effects are a significant issue for RNA-Seq analyses. You can see significant differences in expression due to batch effect.



Adapted from Hicks SC, et al. bioRxiv. 2015. doi:10.1101/025528.

Experimental planning considerations: Batches

How to know whether you have batches?


- Were all RNA isolations performed on the same day?
- Were all library preparations performed on the same day?
- Did the same person perform the RNA isolation/library preparation for all samples?
- Did you use the same reagents for all samples?
- Did you perform the RNA isolation/library preparation in the same location?

If any of the answers was 'No', then you have batches.

Experimental planning considerations: Batches

Best Practice

- Design the experiment in a way to **avoid batches**, if possible.
- If **unable to avoid batches**, **do NOT** confound your experiment by batch




Adapted from Hicks SC, et al. bioRxiv. 2015. doi:10.1101/025528.

Experimental planning considerations: Batches

Best Practice

- If **unable to avoid batches**, **DO** split replicates of the different sample groups across batches
- The more replicates the better (definitely more than 2).



Adapted from Hicks SC, et al. bioRxiv. 2015. doi:10.1101/025528.

Experimental planning considerations: Batches

Best Practice

- If **unable to avoid batches**, **DO** include batch information in your experimental metadata
 - During the analysis, we can **regress out the variation due to batch** so it doesn't affect our results if we have that information.

sample	replicate	condition	batch
sample1	1	control	1
sample2	2	control	1
sample3	3	control	2
sample4	4	control	2
sample5	1	treatment1	1
sample6	2	treatment1	1
sample7	3	treatment1	2
sample8	4	treatment1	2
sample9	1	treatment2	1
sample10	2	treatment2	1
sample11	3	treatment2	2
sample12	4	treatment2	2

These materials have been developed by members of the teaching team at the Harvard Chan Bioinformatics Core (HBC). These are open access materials distributed under the terms of the Creative Commons Attribution license (CC BY 4.0), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited.

